

# 3M™ Surgical Tapes

The full line of latex free and hypoallergenic patient care tapes.



As health care delivery becomes more complex, choosing the right product for an application and using it in the most effective manner is a critical component of successful patient outcomes and cost containment.

Our full range of latex free, hypoallergenic tapes include paper, silk, soft cloth, plastic, foam and cloth adhesive tapes. Each tape has backing and adhesive characteristics for your specific taping technique and situation.

## 3M Surgical Tapes are:

- Latex free
- Hypoallergenic
- Economical

## They are perfect for securing:

- Dressings
- Tubing
- Catheters
- Other devices

**3M Medical**

*from the 3M Health Care family*

**3M**

# 3M family of latex free, hypoallergenic patient care tapes

## 3M™ Transpore™ White Dressing Tape

Gentle tape with strong adhesion and easy, bi-directional tear.

## 3M™ Medipore™ Soft Cloth Surgical Tape

Soft, gentle, conformable tape and pre-cut pad products that are easy to use and easy on the skin.

## 3M™ Medipore™ H Soft Cloth Surgical Tape

The same soft, comfortable tape as standard Medipore tape, but with increased gentleness and stronger adhesion.

## 3M™ Micropore™ Surgical Tape

Gentle, breathable tape with rayon backing. Most trusted brand by clinicians.

## 3M™ Transpore™ Surgical Tape

Transparent, perforated plastic tape that tears easily, works reliably.

## 3M™ Durapore™ Surgical Tape

Silk-like, high strength tape with bi-directional tear.

## 3M™ Microfoam™ Surgical Tape

Highly conformable, elastic foam tape that stretches for compression applications.



## 3M™ Cloth Adhesive Tape

Strong, easy-to-tear tape for securing medical devices and for light immobilization support.

## 3M™ Blenderm™ Surgical Tape

Occlusive, transparent plastic tape that protects from external fluids and contaminants.

## Ordering Information

Cat. No.	Size	Rolls per Box	Boxes per Case
<b>3M™ Transpore™ White Dressing Tape</b>			
1534-0	1/2 in. x 10 yd./1,25 cm x 9,14 m	24	10
1534-1	1 in. x 10 yd./2,5 cm x 9,14 m	12	10
1534-2	2 in. x 10 yd./5,0 cm x 9,14 m	6	10
1534-3	3 in. x 10 yd./7,5 cm x 9,14 m	4	10

## 3M™ Medipore™ Soft Cloth Surgical Tape

2961	1 in. x 10 yd./2,5 cm x 9,14 m	24	1
2962	2 in. x 10 yd./5,0 cm x 9,14 m	12	1
2963	3 in. x 10 yd./7,5 cm x 9,14 m	12	1
2964	4 in. x 10 yd./10 cm x 9,14 m	12	1
2966	6 in. x 10 yd./15 cm x 9,14 m	12	1
2968	8 in. x 10 yd./20 cm x 9,14 m	6	1
2962S	2 in. x 2 yd./5,0 cm x 1,8 m	48	1
2964S	4 in. x 2 yd./10,1 cm x 1,8 m	24	1
2966S	6 in. x 2 yd./15,2 cm x 1,8 m	16	1

## 3M™ Medipore™ H Soft Cloth Surgical Tape

2861	1 in. x 10 yd./2,5 cm x 9,14 m	24	1
2862	2 in. x 10 yd./5,0 cm x 9,14 m	12	1
2863	3 in. x 10 yd./7,5 cm x 9,14 m	12	1
2864	4 in. x 10 yd./10 cm x 9,14 m	12	1
2866	6 in. x 10 yd./15 cm x 9,14 m	12	1
2868	8 in. x 10 yd./20 cm x 9,14 m	6	1
2862S	2 in. x 2 yd./5,0 cm x 1,8 m	48	1
2864S	4 in. x 2 yd./10 cm x 1,8 m	24	1
2866S	6 in. x 2 yd./15 cm x 1,8 m	16	1

## 3M™ Micropore™ Surgical Tape

1530-0	1/2 in. x 10 yd./1,2 cm x 9,14 m	24	10
1530-1	1 in. x 10 yd./2,5 cm x 9,14 m	12	10
1530-2	2 in. x 10 yd./5,0 cm x 9,14 m	6	10
1530-3	3 in. x 10 yd./7,6 cm x 9,14 m	4	10

## Single-use:

1530S-1	1 in. x 1 1/2 yd./2,5 cm x 1,3 m	100	5
1530S-2	2 in. x 1 1/2 yd./5 cm x 1,3 m	50	5

## Tan:

1533-0	1/2 in. 10 yd./1,2 cm x 9,14 m	24	10
1533-1	1 in x 10 yd./2,5 cm 9,14 m	12	10
1533-2	2 in. x 10 yd./5 cm x 9,14 m	6	10

## Dispenser pack:

1535-0	1/2 in. x 10 yd./1,2 cm x 9,14 m	24	10
1535-1	1 in. x 10 yd./2,5 cm x 9,14 m	12	10
1535-2	2 in. x 10 yd./5 cm x 9,14 m	6	10
1535-3	3 in. x 10 yd./7,6 cm x 9,14 m	4	10

Cat. No.	Size	Rolls per Box	Boxes per Case
<b>3M™ Transpore™ Surgical Tape</b>			
1527-0	1/2 in. x 10 yd./1,2 cm x 9,14 m	24	10
1527-1	1 in. x 10 yd./2,5 cm x 9,14 m	12	10
1527-2	2 in. x 10 yd./5 cm x 9,14 m	6	10
1527-3	3 in. x 10 yd./7,6 cm x 9,14 m	4	10

## Single-use:

1527S-1	1 in. x 1 1/2 yd./2,5 cm x 1,3 m	100	5
1527S-2	2 in. x 1 1/2 yd./5 cm x 1,3 m	50	5

## 3M™ Durapore™ Surgical Tape

1538-0	1/2 in. x 10 yd./1,2 cm x 9,14 m	24	10
1538-1	1 in. x 10 yd./2,5 cm x 9,14 m	12	10
1538-2	2 in. x 10 yd./5 cm x 9,14 m	6	10
1538-3	3 in. x 10 yd./7,6 cm 9,14 m	4	10

## Single-use:

1538S-1	1 in. x 1 1/2 yd./2,5 cm x 1,3 m	100	5
1538S-2	2 in. x 1 1/2 yd./5 cm x 1,3 m	50	5

## 3M™ Microfoam™ Surgical Tape

1528-1	1 in. x 5 1/2 yd./2,5 cm x 5 m*	12	6
1528-2	2 in. x 5 1/2 yd./5 cm x 5 m*	6	6
1528-3	3 in. x 5 1/2 yd./7,6 cm x 5 m*	4	6
1528-4	4 in x 5 1/2 yd./10 cm x 5 m*	3	6

## 3M™ Cloth Adhesive Tape

2950-0	1/2 in. x 10 yd./1,25 cm x 9,1 m	24	10
2950-1	1 in. x 10 yd./2,5 cm x 9,1 m	12	10
2950-2	2 in. x 10 yd./5,0 cm x 9,1 m	6	10
2950-3	3 in. x 10 yd./7,6 cm x 9,1 m	4	10

## 3M™ Blenderm™ Surgical Tape

1525-0	1/2 in. x 5 yd./1,25 cm x 4,5 m	24	10
1525-1	1 in. x 5 yd./2,5 cm x 4,5 m	12	10
1525-2	2 in. x 5 yd./5,0 cm x 4,5 m	6	10

\*Stretched length

Cat. No.	Pad Size (3 sheets per pad)	Pads per Box	Boxes per Case
----------	-----------------------------	--------------	----------------

## 3M™ Medipore™ Pre-cut Dressing Covers

2954	3 7/8 in. x 4 7/8 in./10 cm x 12 cm	25	4
2955	3 7/8 in. x 7 7/8 in./10 cm x 20 cm	25	4
2956	5 7/8 in. x 5 7/8 in./15 cm x 15 cm	25	4
2957	5 7/8 in. x 11 in./15 cm x 28 cm	25	4
2958	7 7/8 in. x 11 in./20 cm x 28 cm	25	4

For more information, contact your 3M Health Care sales representative, or call the 3M Health Care Customer Helpline at 1-800-228-3957. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



## Health Care

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA  
1-800-228-3957  
healthcare@3M.com  
www.3M.com/healthcare

## 3M Canada

Post Office Box 5757  
London, Ontario, N6A4T1  
Canada  
1-800-563-2921

3M, Blenderm, Durapore, Medipore, Microfoam, Micropore and Transpore are trademarks of 3M.  
Printed in U.S.A.  
© 3M 2003, Rev 11/06.  
All Rights Reserved.  
70-2009-4185-7 (126.1)ii

Please Recycle