

3M™ Coban™ Self-Adherent Wrap* and 3M™ Coban™ LF Latex Free Self-Adherent Wrap

The self-adherent elastic wrap that functions like a tape,
but sticks only to itself.



1
*Venous Leg
Ulcer Compression*



2
*Surgical Drape
Stabilization*



3
Splint Overwrap



4
*I.V. Catheter
Stabilization*

Features and benefits

- Available non-sterile, sterile, latex free and sterile latex free
- Sticks only to itself; no clips or pins required
- Comfortable for patients
- Lightweight, porous
- Protects primary dressings
- Maintains degree of compression or tension you choose over time
- Available in a variety of colors and sizes (refer to reverse side)

Suggested applications

- Wrapping splints/arm boards
- Compression dressing
- Overwrap for Unna's boot
- Holding other dressings
- Securing I.V.s and other devices especially for combative patients or those with fragile skin
- Sterile drape stabilization and other O.R. applications

3M™ Coban™ Self-Adherent Wrap and 3M™ Coban™ LF Latex Free Self-Adherent Wrap



3M's Skin Health Program provides you with a single source for all your skin health needs including the maintenance of healthy skin, the protection of skin at risk, and the treatment of damaged skin.

To learn more about 3M's Skin Health Program and the Coban family of products, visit us at www.3M.com/skinhealth, contact your 3M Health Care Sales Representative or call the 3M Health Care Customer Helpline at **1-800-228-3957**. Outside of the United States, contact the local 3M subsidiary.



5

Ordering Information

3M™ Coban™ Self-Adherent Wrap*

Catalog Number	Size: Fully-stretched length	Rolls/Case
1581 (Tan)*	1 in. x 5 yd./25 mm x 4,5 m	30
1582 (Tan)*	2 in. x 5 yd./50 mm x 4,5 m	36
1583 (Tan)*	3 in. x 5 yd./75 mm x 4,5 m	24
1584 (Tan)*	4 in. x 5 yd./100 mm x 4,5 m	18
1584L (Tan—Longer Length)*	4 in. x 6.5 yd./100 mm x 6 m	18
1586 (Tan)*	6 in. x 5 yd./150 mm x 4,5 m	12
1583S (Tan—Sterile)*	3 in. x 5 yd./75 mm x 4,5 m	24
1584S (Tan—Sterile)*	4 in. x 5 yd./100 mm x 4,5 m	18
1586S (Tan—Sterile)*	6 in. x 5 yd./150 mm x 4,5 m	12
1581B (Blue)*	1 in. x 5 yd./25 mm x 4,5 m	30
1582B (Blue)*	2 in. x 5 yd./50 mm x 4,5 m	36
1583B (Blue)*	3 in. x 5 yd./75 mm x 4,5 m	24
1583G (Green)*	3 in. x 5 yd./75mm x 4,5 m	24
1583R (Red)*	3 in. x 5 yd./75 mm x 4,5 m	24
1583W (White)*	3 in. x 5 yd./75 mm x 4,5 m	24
1583Y (Yellow)*	3 in. x 5 yd./75 mm x 4,5 m	24
1583A (Color Assortment)*	3 in. x 5 yd./75 mm x 4,5 m	12
Two rolls each: Tan, Blue, Red, Green, Yellow, White		
1583N (Bright Color Assortment)*	3 in. x 5 yd./75 mm x 4,5 m	12
Two rolls each: Bright Green, Bright Orange, Bright Pink; Purple, White, Black		
1584B (Blue)*	4 in. x 5 yd./100 mm x 4,5 m	18

* Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

3M™ Coban™ LF Latex Free Self-Adherent Wrap

Catalog Number	Size: Fully-stretched length	Rolls/Case
2081 (Tan)	1 in. x 5 yd./25 mm x 4,5 m	30
2082 (Tan)	2 in. x 5 yd./50 mm x 4,5 m	36
2083 (Tan)	3 in. x 5 yd./75 mm x 4,5 m	24
2084 (Tan)	4 in. x 5 yd./100 mm x 4,5 m	18
2084L (Tan—Longer Length)	4 in. x 6.5 yd./100 mm x 6 m	18
2086 (Tan)	6 in. x 5 yd./150 mm x 4,5 m	12
2083S (Tan—Sterile)	3 in. x 5 yd./75 mm x 4,5 m	24
2084S (Tan—Sterile)	4 in. x 5 yd./100 mm x 4,5 m	18
2086S (Tan—Sterile)	6 in. x 5 yd./150 mm x 4,5 m	12

Note: Coban wrap and Coban LF latex free wrap are both available in Custom Procedure Trays. Contact your OEM kit manufacturer or 3M Health Care Sales Representative for details.



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1 800 228-3957
healthcare@3M.com
www.3M.com/healthcare

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1 800 563-2921



Minimum 10%
Post-Consumer Fiber

H.I. 4509
Printed in U.S.A.
© 3M 2002 All Rights Reserved
70-2008-5940-6